



DGA Medical, LLC
A Medical Parts & Equipment Company
602-266-6389

A Veteran-Owned Small Business

NEW Biodex Econo Echocardiography Table

Loaded with features...economically priced.

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Allows seated or standing work positions
- Extra wide top with 500-lb patient capacity, suitable for bariatric patients
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Motorized height and Fowler back
- Swing-out cardiac scanning cutout doubles as armrest
- Adjustable sonographer cutout for maximum patient access and comfortable scanning positioning
- Vascular scanning arm board mounts to either side of table

Similar to the Echo Pro with feature and function, the Econo Echocardiography Table is offered without Trendelenburg positions or central locking casters to deliver a lower cost alternative. There are no barriers between sonographer and patient on the Econo Echocardiography Table, allowing correct scanning posture to prevent musculoskeletal injury. The cardiac scanning cutout provides open access to the patient's left thorax area for an unobstructed apical approach. The motorized Fowler back is infinitely adjustable to 80°.

The Econo Echocardiography Table is so accommodating it reduces the total amount of time required to achieve a quality image. Patient comfort is assured even through lengthy procedures by cushioned arm supports, a thick mattress and a unique swing-out cardiac scanning cutout cushion that doubles as an armrest.

Suitable Procedures:

- Abdomen
- Aorta
- Thyroid
- Breast
- Carotid
- Upper Extremity Arterial*
- Renal Transplant
- SMA
- Stress Echo
- Renal
- Pelvis
- Testicular
- Musculoskeletal
- Upper Extremity Venous*
- Renal Artery
- Abdominal Doppler
- Echocardiography
- Transesophageal Echo

* Requires optional vascular arm board



Specifications

- Dimensions: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails
- Height Adjustable: 23" to 39" (58.4 to 99 cm)
- Tabletop: Two primary sections: torso section: features sonographer and patient cutouts with Fowler back positioning, and lower body section remains fixed
 - Cardiac Scanning Cutout: with pivoting cushion that doubles as an armrest
 - Adjustable Sonographer Cutout: drops down or folds up for maximum access to the patient
- Mattress: Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; foam rubber pad covered with durable Naugahyde®
- Upholstery Color: Graphite
- Motions:
 - Height Adjustable: 23" to 39" (58.4 to 99 cm)
 - Fowler Back: 0° to 80° infinitely adjustable
- Controls:
 - Hand Control: Activates height, Fowler back and auto level motions
- Wheels: 5" (12.7 cm) Individual locking swivel casters
- Patient Restraints: Two body straps
- Finish: Powder coat
- Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
- Weight: 335 lb (152 kg)
- Power: 115 VAC or 230 VAC
- Warranty: Two years
- Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN60601-1 standards.



058-701 Table, Echocardiography, Econo, 115 VAC
Includes Hand Controller.
*All options available separately**

058-706 Table, Echocardiography, Econo, 230 VAC
Includes Hand Controller.
*All options available separately**

Optional:

058-632 Side Rails, Folding

058-735 Arm Board, Vascular Scanning

058-656 IV Pole

058- Dispenser, Paper
610 with Cutter Strap

Replacement:

056- Table Paper, Hygienic, 12 rolls
612 18" w x 225' l (46 cm x 69 m)

**Note: Options are installed by the customer.*

US Patent 7,103,932 B1



Biodex™ is a registered trademark of Biodex Medical Systems,
Inc.
© Copyright Biodex Medical Systems, Inc. 2010