



DGA Medical, LLC
A Medical Parts & Equipment Company
602-266-6389

A Veteran-Owned Small Business

NEW Biodex Econo Ultrasound Table

Loaded with features...economically priced.

- Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning
- Allows seated or standing work positions
- 500-lb patient capacity, suitable for bariatric patients
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Motorized height and Fowler back
- Hand controller activates table positions
- Retractable stirrups extend extra long for larger patients
- Optional scanning arm board mounts to either side of table
- Optional side rails fold out of the way to allow unobstructed patient access

Similar to the Ultra Pro with feature and function, the Econo Ultrasound Table is offered without Trendelenburg positions or central locking casters to deliver a lower cost alternative. The Econo Ultrasound Table is a dedicated table for all general ultrasound procedures including OB/GYN. Choose from an assortment of table options to customize the table to suit your department needs, the table is ergonomically designed to suit both patient and sonographer. The table is so accommodating, it reduces the total amount of time required to achieve a quality image.

There are no barriers between sonographer and patient, allowing correct scanning posture to prevent musculoskeletal injury. The open frame design allows sonographers to get closer to the patient, especially useful for carotid imaging or long leg arterial scanning. The motorized Fowler back is infinitely adjustable up to 80° via hand controller. The leg section drops down to 40° for patients with circulatory issues and 80° for stirrup access.

Suitable Procedures:

- Abdomen
- Aorta
- Endovaginal**
- Hysterosonography
- Testicular
- Musculoskeletal
- DVT/Lower Extremity Venous
- Upper Extremity Arterial*
- Renal Transplant
- Renal
- Pelvis
- OB
- Thyroid
- Breast
- Carotid
- Lower Extremity Arterial
- Renal Artery
- Abdominal Doppler



**Requires optional vascular arm board.*

***Requires optional stirrups.*

Specifications

- Dimensions: 70" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with side rails
- Drop-Down Leg Section: 12.6" l (32 cm)
- Height Adjustable: 23" to 39" (58.4 to 99 cm)
- Tabletop: Three primary sections: torso section with Fowler back positioning, center section remains fixed and leg section drops down to 40° and 80° for stirrup access
- Mattress: Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; foam rubber pad covered with durable Naugahyde®
- Upholstery Color: Graphite
- Motions:
 - Height Adjustable: 23" to 39" (58.4 to 99 cm)
 - Fowler Back: 0° to 80° infinitely adjustable
- Controls:
 - Hand Control: Activates height, Fowler back
- Wheels: 5" (12.7 cm) Individual locking swivel casters
- Patient Restraints: Two body straps
- Finish: Powder coat
- Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
- Weight: 335 lb (152 kg)
- Power: 115 VAC or 230 VAC
- Warranty: Two years
- Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN 60601-1 standards.



•

058-726 Table, Ultrasound, Econo, 115 VAC
Includes Hand Controller.
*All options available separately**

058-727 Table, Ultrasound, Econo, 230 VAC
Includes Hand Controller.
*All options available separately**



Optional:



058-632 Side Rails, Folding



058-651 Stirrups, Retractable



058-735 Arm Board, Vascular Scanning




058-656 IV Pole



058-610 Dispenser, Paper
with Cutter Strap

Replacement:



056-612 Table Paper, Hygienic, 12 rolls
18" w x 225' l (46 cm x 69 m)

**Note: Options are installed by the customer.*

