



Accommodates
BARIATRIC
 Patients

Biodex MRI Fowler Back Stretcher

Fowler positioning comfortably accommodates patients for prolonged periods. Ideal for patient transport to and within the MRI suite, the MRI Fowler Back Stretcher is non-magnetic, all aluminum, and designed for heavy-duty use. The Fowler back is infinitely adjustable and the stretcher comes complete with mattress, restraining straps and retractable guard rails. Five-inch ball bearing swivel casters provide easy, controlled mobility and lock firmly to secure the stretcher in position as needed. Optional IV Rod mounts in a socket at the head end of the table.

Fowler back positioning comfortably accommodates patients for prolonged periods

- Fowler positioning
 - Wheels: 5" locking casters
 - Dimensions:
 - Tabletop: 72" l x 24" w (182.8 x 60.9 cm)
 - Overall: 74.5" l x 27.75" w x 34" h (189.3 x 70 x 86 cm)
 - Fowler Back: 0° to 80° infinitely adjustable, manual
 - Patient Capacity: 500 lb (226.8 kg); weight tested to four times the patient load rating.
 - Compliance: MR Conditional; Tested with Philips Achieva 3T System and tested and approved with Siemens Trio 3T System and the Siemens Verio 3T System.
 - Weight: 110 lb (49.9 kg)
 - Warranty: One year
- Finish: Powder coat
 Includes mattress and retractable guard rails
 Shipping Weight: 170 lb (77.1 kg)

240-110 Stretcher, Fowler Back, MRI

240-189 IV Pole 40" H x .70" dia (102 x 1.8cm)

