



DGA MEDICAL LLC
A Medical Equipment & Parts Company
VA FSS # V797P-4278B



Ph 623-251-3599
doug@dgamedical.com
www.dgamedical.com



Biodex MRI Stretcher

This heavy-duty, all aluminum non-magnetic stretcher is an ideal unit for MRI patient transport. Finished in attractive white powder coat, it comes complete with mattress, restraining straps and retractable guard rails. Although light in weight for easy, controlled mobility, the MRI Stretcher accommodates up to a 500-lb patient safely and comfortably.

It rolls on 5" ball-bearing swivel casters that lock firmly to secure the stretcher in place when positioned.

Optional IV Rod mounts in a socket at the head end of the table

Ideal for patient transport, Smooth and effortless mobility, MR conditional for 3T
500-lb patient capacity, Includes mattress and retractable guard rails

- **Dimensions:**
 - Tabletop: 72" l x 24" w (182.8 x 60.9 cm)
 - Overall: 74.5" l x 27.75" w x 34" h (189.3 x 70 x 86 cm)
- **Wheels:** 5" locking casters
- **Finish:** Powder coat
- **Shipping Weight:** 200 lbs.
- **Patient Capacity:** 500 lb (226.8 kg);
weight tested to four times the patient load rating.
- **Compliance:** MR Conditional; Tested with Philips Achieva 3T System and tested and approved with Siemens Trio 3T System and the Siemens Verio 3T System.
- **Weight:** 85 lb (39 kg)
- **Warranty:** One year

240-100 Stretcher, MRI

240-189 IV Pole, 40" H x .70" dia (102 x 1.8cm)

