



DGA Medical, LLC
A Medical Parts & Equipment Company
602-266-6389

Biodex MRI Stretcher

Ideal for patient transport

- Smooth and effortless mobility
- Completely non-ferrous
- 500-lb patient capacity
- Includes mattress and retractable guard rails

This heavy-duty, all aluminum non-magnetic stretcher is an ideal unit for MRI patient transport. Finished in attractive white powder coat, it comes complete with mattress, restraining straps and retractable guard rails. Although light in weight for easy, controlled mobility, the MRI Stretcher accommodates up to a 500-lb patient safely and comfortably.

It rolls on 5" ball-bearing swivel casters that lock firmly to secure the stretcher in place when positioned. Optional IV Rod mounts in a socket at the head end of the table.

Specifications

- Dimensions:
 - Tabletop: 72" l x 24" w (182.8 x 60.9 cm)
 - Overall: 74.5" l x 27.75" w x 34" h (189.3 x 70 x 86 cm)
- Wheels: 5" locking casters
- Finish: Powder coat
- Patient Capacity: 500 lb (226.8 kg); weight tested to four times the patient load rating.
- Compliance: MR Conditional; Tested with Philips Achieva 3T System and tested and approved with Siemens Trio 3T System and the Siemens Verio 3T System.
- Weight: 85 lb (39 kg)
- Shipping Weight: 170 lb (77.1 kg)



240-100 Stretcher, MRI

240-189 IV Pole
40" h x .70" dia (102 x 1.8 cm)

Biodex™ is a registered trademark of Biodex Medical Systems, Inc.
© Copyright Biodex Medical Systems, Inc. 2008

