



DGA MEDICAL LLC
 A Medical Equipment & Parts Company
 VA FSS # V797P-4278B



Ph 623-251-3599
 doug@dgamedical.com
www.dgamedical.com



Biomed Ultra Pro™ Ultrasound Table

The professional's choice when performing general ultrasound procedures, including OB/GYN. Now there is even more reason to choose an Ultra Pro Ultrasound Table.

Enhancements to the ergonomics reduce sonographer fatigue and require less patient repositioning, making it easy to achieve quality images in less time.

The redesigned Fowler back which provides the utmost patient comfort, an enhanced central floor-locking system offering easier access from around the table and an antimicrobial vinyl mattress ensuring safety against infection are all standard on the Ultra Pro.

Customization is made possible with optional items including: flush mounted side rails allowing complete patient access and safety, retractable stirrups with enhanced contouring and extra-long extension for taller patients, an articulating scanning arm board and a headrest extension that drops down 15 degrees to allow for carotid studies.

Suitable Procedures

Abdomen

Hysterosonography

DVT/Lower Extremity Venous*

Renal Transplant

Pelvis

Breast

Upper Extremity Arterial**

Transesophageal Echo

*Requires optional vascular arm board.

**Requires optional stirrups.

Aorta

Testicular

Venous Reflux*

SMA

OB

Carotid

Renal Artery

Endovaginal**

Musculoskeletal

Upper Extremity Venous

Renal

Thyroid

Lower Extremity Arterial

Abdominal Doppler

Ergonomic design for enhanced patient and sonographer comfort

500-lb patient capacity, suitable for bariatric patients

Drop-down, contoured leg section

Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties

Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift

Motorized height, Fowler positioning, Trendelenburg and auto-level adjustment

Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria

Central floor-locking system maintains table stability

Durable hand controller activates table positions

Allows seated or standing work positions

Dimensions:

- Overall: 70" l x 30" w (177.8 x 76.2 cm); 35" w (88.9 cm) with optional side rails
- Drop-Down Leg Section: 12.6" l (32 cm)

Motions:

- Height Adjustable: 23" to 39" (58.4 to 99 cm)
- Trendelenburg: 0° to ±15°
- Fowler Back: 0° to 80° infinitely adjustable

Controls:

- Hand Control: Activates height, Trendelenburg motions,
- Fowler positioning and auto level motions
- Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

Tabletop: Three primary sections: torso section with Fowler positioning, center section remains fixed and leg section drops down to 40° and 80° for stirrup access

Mattress: Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria

Upholstery Color: Graphite

Patient Restraints: Two body straps

Finish: Powder coat

Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.

Weight: 335 lb (152 kg)

Power: 115 VAC or 230 VAC

Warranty: Two-year

Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2

No. 601.1-M90 and EN 60601-1 standards.



Optional Accessories

Foot Controller is ideal for when hands are busy with other functions

Redesigned side rails flush to table for unencumbered patient access

Enhanced stirrups, designed to cup at the heel for added comfort

Articulating Scanning Arm Board, adjustable from 0 to 130 degrees

Headrest extension for easier thyroid and carotid access

IV Pole easily attaches for quick convenience

Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.

All accessories are on VA FSS contract 797P4278B