



Biodex Sound Pro™ Combination Table

Brings together the imaging features of both the Ultra Pro and Echo Prototables. Fully equipped for OB/GYN, general ultrasound and echocardiography procedures, the extra wide tabletop is designed with advanced ergonomics to provide maximum comfort for both patient and sonographer. Designed to conform to natural body extension and sitting position, the motorized Fowler back is infinitely adjustable to 80 degrees via hand or foot control. The table is so accommodating, it reduces the total amount of time required to achieve a quality image. Stress patients are off the treadmill and quickly onto the table in imaging position.

The new cardiac drop-down cushion releases from either side of the table, providing open access to the patient's left thorax area for an unencumbered apical approach. Flush-mounted siderails enable closer access to the patient, without sonographer contortion. When not in use, rails fold beneath the table, eliminating obstruction.

Patient comfort is assured even during lengthy procedures with our new softer feel, Naugahyde® antimicrobial mattress cover, which also provides protection against bacteria. Optional accessories available to expand the function of your table include ergonomic stirrups that cup the heel, an articulating scanning armboard (adjustable from 0 – 130 degrees) and a headrest extension for easy thyroid and carotid access.

Suitable Procedures

Abdomen	Aorta	Endovaginal	Hysterosonography
Testicular	Musculoskeletal	DVT/Lower Extremity Venous	Venous Reflux
Upper Extremity Arterial	Renal Transplant	SMA	Stress Echo
Renal	Pelvis	OB	Thyroid
Breast	Carotid	Lower Extremity Arterial	Upper Extremity
Venous	Renal Artery	Abdominal Doppler	
Echocardiography	Transesophageal Echo		

Features

Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning

Extra wide top with 500-lb patient capacity, suitable for bariatric patients

Cardiac drop-down cushion, allows release and return from either side of table

Industry's largest cardiac scanning cutout – 9.5" w x 15.75" l

Height range down to 23" accommodates wheelchair transfers and hose with ambulatory difficulties

Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift

Motorized height, Fowler positioning, Trendelenburg and auto-level adjustment
Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
Central floor-locking system maintains table stability
Durable hand controller
Allows seated or standing work positions
Redesigned side rails flush to table for unencumbered patient access
Enhanced stirrups, designed to cup at the heel for added comfort
Articulating Scanning Arm Board, adjustable from 0 to 130 degrees• Foot Controller is ideal for when hands are busy with other functions
IV Pole easily attaches for quick convenience
Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.

Dimensions:

- Overall: 71" l x 30" w (180.3 x 76.2 cm); 35" w (88.9 cm) with side rails
- Cardiac Scanning Cutout: 9.5 w x 15.75 l (24.1 x 40 cm)
Drop-down cushion releases from either side of the table.
- Sonographer Cutout: 8" w x 14.75" l (20.3 x 37.4 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.
- Drop-Down Leg Section: 12.6" l (32 cm)

Motions:

- Height Adjustable: 23" to 39" (58.4 to 99 cm)
- Trendelenburg: 0° to ±15°
- Fowler Back: 0° to 80° infinitely adjustable

Controls:

- Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions
 - Foot Control: Activates height, Trendelenburg motions and Fowler positioning
- Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system
Tabletop: Three primary sections: torso section features sonographer and patient cutouts with Fowler positioning; center section remains fixed and leg section drops down to 40° and 80° for stirrup access
- Drop-Down Cardiac Scanning Cushion: Release from either side of the table
 - Adjustable Sonographer Cushion: drops down or folds up for maximum access to the patient

Mattress: Torso section 2" (5 cm) thick, center and leg sections 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria

Upholstery Color: Graphite

Patient Restraints: One body strap

Finish: Powder coat

Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.

Weight: 350 lb (159 kg)

Power: 115 VAC or 230 VAC

Warranty: Two-year

Optional

Headrest extension drops down 15 degrees to allow for carotid and thyroid studies

Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2

No. 601.1-M90 and EN 60601-1 standards.

