



Biodex Echo/Vasc Pro™ Vascular Echocardiography Table

The Echo/Vasc Pro™ Vascular Echocardiography Table is the logical choice for performing vascular scanning and echocardiography. It has been redesigned with many features that help obtain the best quality image and keep the sonographer safe from injury. The new Echo/Vasc Pro has 30 degree reverse Trendelenburg positioning, expanding lower extremity veins for maximum vessel dilation. Enhanced foot support provides the sonographer comfortable access to a patient's calves without crouching. Body straps stabilize the patient and ensure safety throughout the entire procedure. The new design includes a cardiac drop-down cushion, with the ability to release and return from either side of the table, providing ample scanning access for an unobstructed apical approach. As with all Biodex ultrasound tables the Echo/Vasc Pro has enhanced Fowler positioning to accommodate natural body extension and sitting position without slide and shift, redesigned optional side rails for easy access to the patient, an optional Articulating Scanning Arm Board adjustable from 0 to 130-degrees and a new soft antimicrobial vinyl cover.

Suitable Procedures

Echocardiography	Stress Echo	Transesophageal Echo
<u>Lower Extremity</u>		
Venous Ablation	Venous Reflux	Venous Insufficiency
Venous Mapping and Doppler (DVT)		Arterial Scan/Doppler (PAD)

Upper Extremity

Arterial*	Venous Access Placement and Assessment*	
Carotid	Transcranial	Aorta
Breast	Renal Transplant	Testicular
SMA	Carotid Duplex	Abdomen
Thyroid	Renal	Renal Artery

Musculoskeletal

*Requires optional vascular arm board.

All motorized tables are ETL product listed to the UL 60601-1 Medical Electrical Equipment Safety Standard.

Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning

Extra wide top with 500-lb patient capacity, suitable for bariatric patients

Cardiac drop-down cushion, allows release and return from either side of table
Industry's largest cardiac scanning cutout – 9.5" w x 15.75" l
Reverse Trendelenburg to 30 degrees
Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
Motorized height, Fowler positioning, Trendelenburg and auto-level adjustment
Improved foot support with comfortable access to calves without crouching
Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
Central floor-locking system maintains table stability
Durable hand controller
Allows seated or standing work positions
Optional
Foot Controller is ideal for when hands are busy with other functions
Redesigned side rails flush to table for unencumbered patient access
Articulating Scanning Arm Board, adjustable from 0 to 130 degrees
IV Pole easily attaches for quick convenience
Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.
Headrest extension drops down 15 degrees to allow for carotid and thyroid studies
Dimensions:
- Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails
- Cardiac Scanning Cutout: 9.5 w x 15.75 l (24.1 x 40 cm)
Drop-down cushion releases from either side of the table.
- Sonographer Cutout: 8" w x 14.75" l (20.3 x 37.4 cm) Adjustable sonographer cushion drops down or folds up for maximum access to the patient.
Motions:
- Height Adjustable: 23" to 39" (58.4 to 99 cm)
- Trendelenburg: 0° to 15°
- Reverse Trendelenburg: 0° to 30°
- Fowler Back: 0° to 80° infinitely adjustable
Controls:
- Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions
- Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning
Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system
Tabletop: Two primary sections: torso section features sonographer and patient cutouts with Fowler positioning, and lower body section remains fixed
Mattress: Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
Upholstery Color: Graphite
Patient Restraints: Two body straps
Finish: Powder coat
Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
Weight: 335 lb (152 kg)
Power: 115 VAC or 230 VAC
Warranty: Two-year
Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2 No. 601.1-M90 and EN60601-1 standards.
058-702 Table, Vascular Echocardiography, Echo/Vasc Pro™, 115 VAC

Includes Hand Controller and Foot Support.

Optional:

058-740 Foot Controller

058-633 Side Rails, Folding

058-736 Arm Board, Articulating Scanning

058-737 IV Pole

058-611 Dispenser, Paper with Cutter Strap

058-738 Extension, Headrest