



Biodex Vasc Pro™ Vascular Ultrasound Table

The new Vasc Pro™ Vascular Ultrasound Table is designed with ergonomics in mind. It addresses sonographer comfort and safety and proper patient positioning, virtually eliminating repetitive-strain injuries from vascular scanning.

Sonographer comfort and safety are accomplished through new features such as the enhanced foot support allowing easier access to the patient calves without crouching and optional side rails that enable closer access to the patient without sonographer contortion. Fowler positioning to 80 degrees, which accommodates natural body extension and sitting position without slide or shift, an optional Articulating Scanning Arm Board adjustable from 0 to 130-degrees locking every 10 degrees, the thick antimicrobial mattress and table adjustability all create a comfortable imaging experience for the patient.

Suitable Procedures

Lower Extremity

Venous Reflux	Venous Insufficiency	Venous Ablation
Venous Mapping and Doppler (DVT)		Arterial Scan/Doppler (PAD)

Upper Extremity

Arterial*	Carotid	Transcranial
Aorta	Breast	Renal Transplant
Testicular	SMA	Carotid Duplex
Abdomen	Thyroid	Renal
Renal Artery	Musculoskeletal	Venous Access Placement and Assessment*

*Requires optional vascular arm board.

All motorized tables are ETL product listed to the UL 60601-1 Medical Electrical Equipment Safety Standard.

Ergonomic design helps prevent musculoskeletal injuries – comfortable shoulder, arm and hand positions when scanning

Extra wide top with 500-lb patient capacity, suitable for bariatric patients

Reverse Trendelenburg to 30 degrees

Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties

Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift

Motorized height, Fowler positioning, Trendelenburg and auto-level adjustment

Improved foot support with comfortable access to calves without crouching

Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria

Central floor-locking system maintains table stability

Durable hand controller

Allows seated or standing work positions

Optional

Foot Controller is ideal for when hands are busy with other functions

Redesigned side rails flush to table for unencumbered patient access

Articulating Scanning Arm Board, adjustable from 0 to 130 degrees

IV Pole easily attaches for quick convenience

Paper Dispenser accommodates a roll of hygienic table paper. Includes a cutter strap.

Headrest extension drops down 15 degrees to allow for carotid and thyroid studies

Dimensions:

- Overall: 73" l x 30" w (185.4 x 76.2 cm); 35" w (88.9 cm) with optional side rails

Motions:

- Height Adjustable: 23" to 39" (58.4 to 99 cm)

- Trendelenburg: 0° to 15°

- Reverse Trendelenburg: 0° to 30°

- Fowler Back: 0° to 80° infinitely adjustable

Controls:

- Hand Control: Activates height, Trendelenburg motions, Fowler positioning and auto level motions

- Foot Control (optional): Activates height, Trendelenburg motions, and Fowler positioning

Wheels: 5" (12.7 cm) Individual locking swivel casters; central floor-locking system

Tabletop: Two primary sections: torso section with Fowler positioning; lower body section remains fixed

Mattress: Torso section 2" (5 cm) thick, lower body section 3" (8 cm) thick; Naugahyde®, antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria

Upholstery Color: Graphite

Patient Restraints: Three body straps

Finish: Powder coat

Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.

Weight: 335 lb (152 kg)

Power: 115 VAC or 230 VAC

Warranty: Two-year

Certifications: ETL and cETL Listed to UL60601-1 and CAN/CSA C22.2

No. 601.1-M90 and EN60601-1 standards.

058-732 Table, Vascular Ultrasound, Vasc Pro™, 115 VAC,

Includes Hand Controller and Foot Support.

Optional:

058-740 Foot Controller

058-633 Side Rails, Folding

058-736 Arm Board, Articulating Scanning

058-737 IV Pole

058-611 Dispenser, Paper with Cutter Strap

058-738 Extension, Headrest