



Biodex MRI Fowler Back Stretcher

Fowler positioning comfortably accommodates patients for prolonged periods. Ideal for patient transport to and within the MRI suite, the MRI Fowler Back Stretcher is non-magnetic, all aluminum, and designed for heavy-duty use. The Fowler back is infinitely adjustable and the stretcher comes complete with mattress, restraining straps and retractable guard rails. Five-inch ball bearing swivel casters provide easy, controlled mobility and lock firmly to secure the stretcher in position as needed. Optional IV Rod mounts in a socket at the head end of the table.

Fowler back positioning comfortably accommodates patients for prolonged periods

- Fowler positioning
 - Wheels: 5" locking casters
 - Dimensions:
 - Tabletop: 72" l x 24" w (182.8 x 60.9 cm)
 - Overall: 74.5" l x 27.75" w x 34" h (189.3 x 70 x 86 cm)
 - Fowler Back: 0° to 80° infinitely adjustable, manual
 - Patient Capacity: 500 lb (226.8 kg); weight tested to four times the patient load rating.
 - Compliance: MR Conditional; Tested with Philips Achieva 3T System and tested and approved with Siemens Trio 3T System and the Siemens Verio 3T System.
 - Weight: 110 lb (49.9 kg)
 - Warranty: One year
- Finish: Powder coat
Includes mattress and retractable guard rails
Shipping Weight: 170 lb (77.1 kg)

240-110 Stretcher, Fowler Back, MRI

240-189 IV Pole 40" H x .70" dia (102 x 1.8cm)

